



SPXN Polyclonal Antibody

Cat No: HR1AP12532

For research use only

Overview

Product Name	SPXN Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	12kDa
GeneID?Human?	C12orf39
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	spexin hormone(SPX) Homo sapiens The protein encoded by this gene is a hormone involved in modulation of cardiovascular and renal function. It has also been shown in rats to cause weight loss. Several transcript variants have been found for this gene. [provided by RefSeq, Feb 2016],