

# IRGM Polyclonal Antibody

Cat No: HR1AP12525

For research use only

## Overview

Product Name	IRGM Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	19kDa
GeneID?Human?	IRGM IFI1 IRGM1 LRG47
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	immunity related GTPase M(IRGM) Homo sapiens This gene encodes a member of the p47 immunity-related GTPase family. The encoded protein may play a role in the innate immune response by regulating autophagy formation in response to intracellular pathogens. Polymorphisms that affect the normal expression of this gene are associated with a susceptibility to Crohn's disease and tuberculosis.[provided by RefSeq, Oct 2010],