

# RTN2 Polyclonal Antibody

Cat No: HR1AP12522

For research use only

## Overview

Product Name	RTN2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	59kDa
GeneID?Human?	RTN2 NSPL1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	reticulon 2(RTN2) Homo sapiens This gene belongs to the family of reticulon encoding genes. Reticulons are necessary for proper generation of tubular endoplasmic reticulum and likely play a role in intracellular vesicular transport. Alternatively spliced transcript variants encoding different isoforms have been identified. Mutations at this locus have been associated with autosomal dominant spastic paraplegia-12. [provided by RefSeq, Apr 2012],