

TBX20 Polyclonal Antibody

Cat No: HR1AP12519

For research use only

Overview

Product Name	TBX20 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	49kDa
GeneID?Human?	TBX20
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>T-box 20(TBX20) Homo sapiens This gene encodes a T-box family member. The T-box family members share a common DNA binding domain, termed the T-box, and they are transcription factors involved in the regulation of developmental processes. This gene is essential for heart development. Mutations in this gene are associated with diverse cardiac pathologies, including defects in septation, valvulogenesis and cardiomyopathy. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2009].</p>