

FLVC1 Polyclonal Antibody

Cat No: HR1AP12515

For research use only

Overview

Product Name	FLVC1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	61kDa
GeneID?Human?	FLVCR1 FLVCR
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	feline leukemia virus subgroup C cellular receptor 1 (FLVCR1) Homo sapiens This gene encodes a member of the major facilitator superfamily of transporter proteins. The encoded protein is a heme transporter that may play a critical role in erythropoiesis by protecting developing erythroid cells from heme toxicity. This gene may play a role in posterior column ataxia with retinitis pigmentosa and the hematological disorder Diamond-Blackfan syndrome. [provided by RefSeq, Jan 2011],