

TGFA1 Polyclonal Antibody

Cat No: HR1AP12513

For research use only

Overview

Product Name	TGFA1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	94kDa
GeneID?Human?	TGFBRAP1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	transforming growth factor beta receptor associated protein 1(TGFBRAP1) Homo sapiens This gene encodes a protein that binds to transforming growth factor-beta (TGF-beta) receptors and plays a role in TGF-beta signaling. The encoded protein acts as a chaperone in signaling downstream of TGF-beta. It is involved in signal-dependent association with SMAD4. The protein is also a component of mammalian CORVET, a multisubunit tethering protein complex that is involved in fusion of early endosomes. [provided by RefSeq, Jun 2016].