

PF4V Polyclonal Antibody

Cat No: HR1AP12510

For research use only

Overview

Product Name	PF4V Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	11kDa
GeneID?Human?	PF4V1 CXCL4V1 SCYB4V1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	platelet factor 4 variant 1(PF4V1) Homo sapiens The protein encoded by this gene is a chemokine that is highly similar to platelet factor 4. The encoded protein displays a strong antiangiogenic function and is regulated by chemokine (C-X-C motif) receptor 3. This protein also impairs tumor growth and can protect against blood-retinal barrier breakdown in diabetes patients. [provided by RefSeq, Nov 2015],