

# SCTM1 Polyclonal Antibody

Cat No: HR1AP12511

For research use only

## Overview

Product Name	SCTM1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	27kDa
GeneID?Human?	SECTM1 K12
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	secreted and transmembrane 1(SECTM1) Homo sapiens This gene encodes a transmembrane and secreted protein with characteristics of a type 1a transmembrane protein. It is found in a perinuclear Golgi-like pattern and thought to be involved in hematopoietic and/or immune system processes. [provided by RefSeq, Jul 2008],