

TRIP6 Polyclonal Antibody

Cat No: HR1AP12512

For research use only

Overview

Product Name	TRIP6 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	52kDa
GeneID?Human?	TRIP6 OIP1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>thyroid hormone receptor interactor 6(TRIP6) Homo sapiens This gene is a member of the zyxin family and encodes a protein with three LIM zinc-binding domains. This protein localizes to focal adhesion sites and along actin stress fibers. Recruitment of this protein to the plasma membrane occurs in a lysophosphatidic acid (LPA)-dependent manner and it regulates LPA-induced cell migration. Alternatively spliced variants which encode different protein isoforms have been described; however, not all variants have been fully characterized. [provided by RefSeq, Jul 2008],</p>