

ID3 Polyclonal Antibody

Cat No: HR1AP12509

For research use only

Overview

Product Name	ID3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	13kDa
GeneID?Human?	ID3 1R21 BHLHB25 HEIR1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	inhibitor of DNA binding 3, HLH protein(ID3) Homo sapiens The protein encoded by this gene is a helix-loop-helix (HLH) protein that can form heterodimers with other HLH proteins. However, the encoded protein lacks a basic DNA-binding domain and therefore inhibits the DNA binding of any HLH protein with which it interacts. [provided by RefSeq, Aug 2011],