

## UPK3A Polyclonal Antibody

Cat No: HR1AP12507

For research use only

### Overview

Product Name	UPK3A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	31kDa
GeneID?Human?	UPK3A UPK3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>uroplakin 3A(UPK3A) Homo sapiens This gene encodes a member of the uroplakin family, a group of transmembrane proteins that form complexes on the apical surface of the bladder epithelium. Mutations in this gene may be associated with renal adysplasia. Alternatively spliced transcript variants have been described.[provided by RefSeq, Nov 2009],</p>