

SALL4 Polyclonal Antibody

Cat No: HR1AP12505

For research use only

Overview

Product Name	SALL4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	115kDa
GeneID?Human?	SALL4 ZNF797
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	spalt like transcription factor 4(SALL4) Homo sapiens This gene encodes a zinc finger transcription factor thought to play a role in the development of abducens motor neurons. Defects in this gene are a cause of Duane-radial ray syndrome (DRRS). Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Dec 2015],