

# ES1 Polyclonal Antibody

Cat No: HR1AP12503

For research use only

## Overview

Product Name	ES1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	29kDa
GeneID?Human?	C21orf33 HES1 KNPI
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	chromosome 21 open reading frame 33(C21orf33) Homo sapiens This gene encodes a potential mitochondrial protein that is a member of the DJ-1/Pfpl gene family. This protein is overexpressed in fetal Down syndrome brain. Alternate splicing results in multiple transcript variants. [provided by RefSeq, May 2010],