



NRG4 Polyclonal Antibody

Cat No: HR1AP12497

For research use only

Overview

Product Name	NRG4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	12kDa
GeneID?Human?	NRG4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	neuregulin 4(NRG4) Homo sapiens The neuregulins, including NRG4, activate type-1 growth factor receptors (see EGFR; MIM 131550) to initiating cell-to-cell signaling through tyrosine phosphorylation (Harari et al., 1999 [PubMed 10348342]).[supplied by OMIM, Mar 2008],