

SLUR1 Polyclonal Antibody

Cat No: HR1AP12496

For research use only

Overview

| | |
|-----------------------|---|
| Product Name | SLUR1 Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,ELISA |
| Species Reactivity | Human |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 11kDa |
| GeneID?Human? | SLURP1 ARS |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | |
| Background | secreted LY6/PLAUR domain containing 1(SLURP1) Homo sapiens The protein encoded by this gene is a member of the Ly6/uPAR family but lacks a GPI-anchoring signal sequence. It is thought that this secreted protein contains antitumor activity. Mutations in this gene have been associated with Mal de Meleda, a rare autosomal recessive skin disorder. This gene maps to the same chromosomal region as several members of the Ly6/uPAR family of glycoprotein receptors. [provided by RefSeq, Jul 2008], |