

TDGF1 Polyclonal Antibody

Cat No: HR1AP12491

For research use only

Overview

Product Name	TDGF1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	20kDa
GeneID?Human?	TDGF1 CRIPTO
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>teratocarcinoma-derived growth factor 1(TDGF1) Homo sapiens This gene encodes an epidermal growth factor-related protein that contains a cripto, FRL-1, and cryptic domain. The encoded protein is an extracellular, membrane-bound signaling protein that plays an essential role in embryonic development and tumor growth. Mutations in this gene are associated with forebrain defects. Pseudogenes of this gene are found on chromosomes 2, 3, 6, 8, 19 and X. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Mar 2010],</p>