

RNAS4 Polyclonal Antibody

Cat No: HR1AP12478

For research use only

Overview

Product Name	RNAS4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	16kDa
GeneID?Human?	RNASE4 RNS4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>ribonuclease A family member 4(RNASE4) Homo sapiens The protein encoded by this gene belongs to the pancreatic ribonuclease family. It plays an important role in mRNA cleavage and has marked specificity towards the 3' side of uridine nucleotides. Alternative splicing results in four transcript variants encoding the same protein. This gene and the gene that encodes angiogenin share promoters and 5' exons. Each gene splices to a unique downstream exon that contains its complete coding region. [provided by RefSeq, Aug 2013],</p>