

YIPF3 Polyclonal Antibody

Cat No: HR1AP12479

For research use only

Overview

Product Name	YIPF3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	38kDa
GeneID?Human?	YIPF3 C6orf109 KLIP1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	developmental stage:Expressed in fetal liver (at protein level). Expressed in embryonic hematopoiesis sites.,function:May play a role in hematopoiesis.,PTM:Not glycosylated.,similarity:Belongs to the YIP1 family.,subcellular location:Localization to the cytoplasm or to the cell membrane is developmentally and ontogenetically regulated.,tissue specificity:Expressed by nucleated hematopoietic cells (at protein level).