

CREL1 Polyclonal Antibody

Cat No: HR1AP12476

For research use only

Overview

Product Name	CREL1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	46kDa
GeneID?Human?	CRELD1 CIRRIN UNQ188/PRO214
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	cysteine rich with EGF like domains 1(CRELD1) Homo sapiens This gene encodes a member of a subfamily of epidermal growth factor-related proteins. The encoded protein is characterized by a cysteine-rich with epidermal growth factor-like domain. This protein may function as a cell adhesion molecule. Mutations in this gene are the cause of atrioventricular septal defect. Alternate splicing results in multiple transcript variants.[provided by RefSeq, Apr 2010].