

CHODL Polyclonal Antibody

Cat No: HR1AP12473

For research use only

Overview

Product Name	CHODL Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	30kDa
GeneID?Human?	CHODL C21orf68 PRED12 UNQ872/PRO1890
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>chondrolectin(CHODL) Homo sapiens This gene encodes a type I membrane protein with a carbohydrate recognition domain characteristic of C-type lectins in its extracellular portion. In other proteins, this domain is involved in endocytosis of glycoproteins and exogenous sugar-bearing pathogens. This protein localizes predominantly to the perinuclear region. Several transcript variants encoding a few different isoforms have been found for this gene. [provided by RefSeq, Feb 2011].</p>