

MEP1A Polyclonal Antibody

Cat No: HR1AP12469

For research use only

Overview

Product Name	MEP1A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	82kDa
GeneID?Human?	MEP1A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>catalytic activity:Hydrolysis of protein and peptide substrates preferentially on carboxyl side of hydrophobic residues.,cofactor:Binds 1 zinc ion per subunit.,similarity:Belongs to the peptidase M12A family.,similarity:Contains 1 EGF-like domain.,similarity:Contains 1 MAM domain.,similarity:Contains 1 MATH domain.,subunit:Homotetramer of alpha or beta subunits; heterotetramer of two alpha and two beta subunits are formed by non-covalent association of two disulfide-linked heterodimers.,</p>