

SMC1B Polyclonal Antibody

Cat No: HR1AP12463

For research use only

Overview

Product Name	SMC1B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	135kDa
GeneID?Human?	SMC1B SMC1L2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	structural maintenance of chromosomes 1B(SMC1B) Homo sapiens SMC1L2 belongs to a family of proteins required for chromatid cohesion and DNA recombination during meiosis and mitosis (3:Revenkova et al., 2001 [PubMed 11564881]).[supplied by OMIM, Mar 2008],