

BACE2 Polyclonal Antibody

Cat No: HR1AP12459

For research use only

Overview

Product Name	BACE2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	56kDa
GeneID?Human?	BACE2 AEPLC ALP56 ASP21 CDA13 UNQ418/PRO852
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>beta-site APP-cleaving enzyme 2(BACE2) Homo sapiens This gene encodes an integral membrane glycoprotein that functions as an aspartic protease. The encoded protein cleaves amyloid precursor protein into amyloid beta peptide, which is a critical step in the etiology of Alzheimer's disease and Down syndrome. The protein precursor is further processed into an active mature peptide. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2013].</p>