

# AXDN1 Polyclonal Antibody

Cat No: HR1AP12457

For research use only

## Overview

Product Name	AXDN1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	111kDa
GeneID?Human?	AXDND1 C1orf125
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	