

RTN4R Polyclonal Antibody

Cat No: HR1AP12458

For research use only

Overview

Product Name	RTN4R Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	52kDa
GeneID?Human?	RTN4R NOGOR UNQ330/PRO526
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	reticulon 4 receptor(RTN4R) Homo sapiens This gene encodes the receptor for reticulon 4, oligodendrocyte myelin glycoprotein and myelin-associated glycoprotein. This receptor mediates axonal growth inhibition and may play a role in regulating axonal regeneration and plasticity in the adult central nervous system. [provided by RefSeq, Jul 2008],