

PARPT Polyclonal Antibody

Cat No: HR1AP12452

For research use only

Overview

Product Name	PARPT Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	72kDa
GeneID?Human?	TIPARP PARP7
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	TCDD inducible poly(ADP-ribose) polymerase(TIPARP) Homo sapiens This gene encodes a member of the poly(ADP-ribose) polymerase superfamily. Studies of the mouse ortholog have shown that the encoded protein catalyzes histone poly(ADP-ribosyl)ation and may be involved in T-cell function. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2010],