

## CC14C Polyclonal Antibody

Cat No: HR1AP12448

For research use only

### Overview

Product Name	CC14C Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	60kDa
GeneID?Human?	CDC14C CDC14B2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	Dual-specificity phosphatase. Preferentially dephosphorylates proteins modified by proline-directed kinases (By similarity).