

F19A5 Polyclonal Antibody

Cat No: HR1AP12446

For research use only

Overview

Product Name	F19A5 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	14kDa
GeneID?Human?	FAM19A5 TAF5 UNQ5208/PRO34524
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>family with sequence similarity 19 member A5, C-C motif chemokine like(FAM19A5) Homo sapiens This gene is a member of the TAF5 family which is composed of five highly homologous genes that encode small secreted proteins. These proteins contain conserved cysteine residues at fixed positions, and are distantly related to MIP-1alpha, a member of the CC-chemokine family. The TAF5 proteins are predominantly expressed in specific regions of the brain, and are postulated to function as brain-specific chemokines or neurokinins that act as regulators of immune and nervous cells. [provided by RefSeq, Sep 2013],</p>