



SNAI3 Polyclonal Antibody

Cat No: HR1AP12441

For research use only

Overview

Product Name	SNAI3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	32kDa
GeneID?Human?	SNAI3 ZNF293
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	snail family transcriptional repressor 3(SNAI3) Homo sapiens SNAI3 is a member of the SNAIL gene family, named for the Drosophila snail gene, which plays roles in mesodermal formation during embryogenesis (Katoh and Katoh, 2003 [PubMed 12579345]).[supplied by OMIM, Apr 2009],