

# MFAP5 Polyclonal Antibody

Cat No: HR1AP12439

For research use only

## Overview

Product Name	MFAP5 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	19kDa
GeneID?Human?	MFAP5 MAGP2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	microfibrillar associated protein 5(MFAP5) Homo sapiens This gene encodes a 25-kD microfibril-associated glycoprotein which is a component of microfibrils of the extracellular matrix. The encoded protein promotes attachment of cells to microfibrils via alpha-V-beta-3 integrin. Deficiency of this gene in mice results in neutropenia. Alternate splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2014],