

## FCRLB Polyclonal Antibody

Cat No: HR1AP12436

For research use only

### Overview

Product Name	FCRLB Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	46kDa
GeneID?Human?	FCRLB FCRL2 FCRLM2 FCRY FREB2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	Fc receptor like B(FCRLB) Homo sapiens FCRL2 belongs to the Fc receptor family. Fc receptors are involved in phagocytosis, antibody-dependent cell cytotoxicity, immediate hypersensitivity, and transcytosis of immunoglobulins via their ability to bind immunoglobulin (Ig) constant regions (Chikaev et al., 2005 [PubMed 15676285]).[supplied by OMIM, Mar 2008],