

## HECAM Polyclonal Antibody

Cat No: HR1AP12434

For research use only

### Overview

Product Name	HECAM Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	45kDa
GeneID?Human?	HEPACAM
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>hepatic and glial cell adhesion molecule(HEPACAM) Homo sapiens The protein encoded by this gene is a single-pass type I membrane protein that localizes to the cytoplasmic side of the cell membrane. The encoded protein acts as a homodimer and is involved in cell motility and cell-matrix interactions. The expression of this gene is downregulated or undetectable in many cancer cell lines, so this may be a tumor suppressor gene. [provided by RefSeq, Jul 2011],</p>