

## GLMN Polyclonal Antibody

Cat No: HR1AP12430

For research use only

### Overview

Product Name	GLMN Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	65kDa
GeneID?Human?	GLMN FAP48 FAP68 VMGLOM
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>glomulin, FKBP associated protein(GLMN) Homo sapiens This gene encodes a phosphorylated protein that is a member of a Skp1-Cullin-F-box-like complex. The protein is essential for normal development of the vasculature and mutations in this gene have been associated with glomuvenous malformations, also called glomangiomas. Multiple splice variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2016],</p>