

CNDP1 Polyclonal Antibody

Cat No: HR1AP12421

For research use only

Overview

Product Name	CNDP1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	55kDa
GeneID?Human?	CNDP1 CN1 CPGL2 UNQ1915/PRO4380
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>carnosine dipeptidase 1(CNDP1) Homo sapiens This gene encodes a member of the M20 metalloprotease family. The encoded protein is specifically expressed in the brain, is a homodimeric dipeptidase which was identified as human carnosinase. This gene contains trinucleotide (CTG) repeat length polymorphism in the coding region. [provided by RefSeq, Jul 2008],</p>