

EGFL7 Polyclonal Antibody

Cat No: HR1AP12416

For research use only

Overview

Product Name	EGFL7 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	30kDa
GeneID?Human?	EGFL7 MEGF7 UNQ187/PRO1449
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	EGF like domain multiple 7(EGFL7) Homo sapiens This gene encodes a secreted endothelial cell protein that contains two epidermal growth factor-like domains. The encoded protein may play a role in regulating vasculogenesis. This protein may be involved in the growth and proliferation of tumor cells. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Feb 2012],