

WNT2 Polyclonal Antibody

Cat No: HR1AP12411

For research use only

Overview

Product Name	WNT2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	39kDa
GeneID?Human?	WNT2 INT1L1 IRP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>Wnt family member 2(WNT2) Homo sapiens This gene is a member of the WNT gene family. The WNT gene family consists of structurally related genes which encode secreted signaling proteins. These proteins have been implicated in oncogenesis and in several developmental processes, including regulation of cell fate and patterning during embryogenesis. Alternatively spliced transcript variants have been identified for this gene. [provided by RefSeq, Jul 2008],</p>