

# SLAF8 Polyclonal Antibody

Cat No: HR1AP12408

For research use only

## Overview

Product Name	SLAF8 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	31kDa
GeneID?Human?	SLAMF8 BLAME
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	SLAM family member 8(SLAMF8) Homo sapiens This gene encodes a member of the CD2 family of cell surface proteins involved in lymphocyte activation. These proteins are characterized by Ig domains. This protein is expressed in lymphoid tissues, and studies of a similar protein in mouse suggest that it may function during B cell lineage commitment. The gene is found in a region of chromosome 1 containing many CD2 genes. [provided by RefSeq, Jul 2008],