

R4RL1 Polyclonal Antibody

Cat No: HR1AP12404

For research use only

Overview

Product Name	R4RL1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	48kDa
GeneID?Human?	RTN4RL1 NGRH2 NGRL2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:May play a role in regulating axonal regeneration and plasticity in the adult central nervous system.,similarity:Belongs to the Nogo receptor family.,similarity:Contains 8 LRR (leucine-rich) repeats.,subcellular location:Localized to the surface of neurons, including axons.,tissue specificity:Predominantly expressed in brain. Expressed at lower levels in kidney, lung, mammary gland, placenta, salivary gland, skeletal muscle and spleen.,