

# SLIK1 Polyclonal Antibody

Cat No: HR1AP12402

For research use only

## Overview

Product Name	SLIK1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	76kDa
GeneID?Human?	SLITRK1 KIAA1910 LRRC12 UNQ233/PRO266
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>SLIT and NTRK like family member 1(SLITRK1) Homo sapiens This gene encodes a member of the SLITRK protein family. Members of this family are integral membrane proteins that are characterized by two N-terminal leucine-rich repeat (LRR) domains and a C-terminal region that shares homology with trk neurotrophin receptors. However, the protein encoded by this gene lacks the region of homology to neurotrophin receptors. This protein is thought to be involved in neurite outgrowth. Mutations in this gene may be associated with Tourette syndrome. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2013],</p>