



CRIS2 Polyclonal Antibody

Cat No: HR1AP12397

For research use only

Overview

Product Name	CRIS2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	26kDa
GeneID?Human?	CRISP2 GAPDL5 TPX1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:May regulate some ion channels' activity and thereby regulate calcium fluxes during sperm capacitation.,similarity:Belongs to the CRISP family.,tissue specificity:Testis and epididymis.,