

HPLN4 Polyclonal Antibody

Cat No: HR1AP12399

For research use only

Overview

Product Name	HPLN4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	44kDa
GeneID?Human?	HAPLN4 BRAL2 KIAA1926
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function: Binds to hyaluronic acid and may be involved in formation of the extracellular matrix., similarity: Belongs to the HAPLN family., similarity: Contains 1 Ig-like C2-type (immunoglobulin-like) domain., similarity: Contains 2 Link domains., tissue specificity: Expressed predominantly in brain.,