

# LPIN2 Polyclonal Antibody

Cat No: HR1AP12395

For research use only

## Overview

Product Name	LPIN2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	98kDa
GeneID?Human?	LPIN2 KIAA0249
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	lipin 2(LPIN2) Homo sapiens Mouse studies suggest that this gene functions during normal adipose tissue development and may play a role in human triglyceride metabolism. This gene represents a candidate gene for human lipodystrophy, characterized by loss of body fat, fatty liver, hypertriglyceridemia, and insulin resistance. [provided by RefSeq, Jul 2008],