

SLIK4 Polyclonal Antibody

Cat No: HR1AP12393

For research use only

Overview

Product Name	SLIK4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	92kDa
GeneID?Human?	SLITRK4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>SLIT and NTRK like family member 4(SLITRK4) Homo sapiens This gene encodes a transmembrane protein belonging to the the SLITRK family. These family members include two N-terminal leucine-rich repeat domains similar to those found in the axonal growth-controlling protein SLIT, as well as C-terminal regions similar to neurotrophin receptors. Studies of an homologous protein in mouse suggest that this family member functions to suppress neurite outgrowth. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2010].</p>