

## SEM6A Polyclonal Antibody

Cat No: HR1AP12389

For research use only

### Overview

Product Name	SEM6A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	113kDa
GeneID?Human?	SEMA6A KIAA1368 SEMAQ
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	semaphorin 6A(SEMA6A) Homo sapiens The transmembrane semaphorin SEMA6A is expressed in developing neural tissue and is required for proper development of the thalamocortical projection (Leighton et al., 2001 [PubMed 11242070]).[supplied by OMIM, Feb 2011],