

# TRML2 Polyclonal Antibody

Cat No: HR1AP12390

For research use only

## Overview

Product Name	TRML2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	35kDa
GeneID?Human?	TREML2 C6orf76 TLT2 UNQ6268/PRO20473
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	triggering receptor expressed on myeloid cells like 2(TREML2) Homo sapiens TREML2 is located in a gene cluster on chromosome 6 with the single Ig variable (IgV) domain activating receptors TREM1 (MIM 605085) and TREM2 (MIM 605086), but it has distinct structural and functional properties (Allcock et al., 2003 [PubMed 12645956]).[supplied by OMIM, Mar 2008],