

PTHD2 Polyclonal Antibody

Cat No: HR1AP12386

For research use only

Overview

Product Name	PTHD2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	153kDa
GeneID?Human?	PTCHD2 DISP3 KIAA1337
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	developmental stage:Expressed in fetal brain.,similarity:Belongs to the patched family.,similarity:Contains 1 SSD (sterol-sensing) domain.,tissue specificity:Expressed in brain and testis.,