

# IGS11 Polyclonal Antibody

Cat No: HR1AP12384

For research use only

## Overview

Product Name	IGS11 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	47kDa
GeneID?Human?	IGSF11 BTIGSF CXADRL1 VSIG3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	immunoglobulin superfamily member 11(IGSF11) Homo sapiens IGSF11 is an immunoglobulin (Ig) superfamily member that is preferentially expressed in brain and testis. It shares significant homology with coxsackievirus and adenovirus receptor (CXADR; MIM 602621) and endothelial cell-selective adhesion molecule (ESAM).[supplied by OMIM, Apr 2005],