

MYCL1 Polyclonal Antibody

Cat No: HR1AP12380

For research use only

Overview

Product Name	MYCL1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	40kDa
GeneID?Human?	MYCL1 BHLHE38 LMYC MYCL
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Contains 1 basic helix-loop-helix (bHLH) domain.,subunit:Efficient DNA binding requires dimerization with another bHLH protein. Binds DNA as an heterodimer with MAX.,