

NANO2 Polyclonal Antibody

Cat No: HR1AP12381

For research use only

Overview

Product Name	NANO2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	15kDa
GeneID?Human?	NANOS2 NOS2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>domain:The Nanos-type zinc finger is composed of two C2HC motifs, each motif binding one molecule of zinc. It is essential for the translation repression activity of the protein.,function:Required to support proliferation and self-renewal of proximal germ cells in males only. Probably regulates translation of specific mRNAs by associating with the 3'-UTR of mRNA targets. Essential for spermatogonia formation.,similarity:Belongs to the nanos family.,similarity:Contains 1 nanos-type zinc finger.,</p>