

CDX4 Polyclonal Antibody

Cat No: HR1AP12378

For research use only

Overview

Product Name	CDX4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	31kDa
GeneID?Human?	CDX4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	caudal type homeobox 4(CDX4) Homo sapiens This gene encodes a member of a small subfamily of homeobox containing transcription factors. The encoded protein may regulate homeobox gene expression during anteroposterior patterning and hematopoiesis. [provided by RefSeq, Aug 2012],